

The Methacholine Challenge Test for Reversible Airways Disease Assessment: A Practical Guide on how to Interpret New 2017 ERS Guidelines

Jason A Suggett¹, Mark W Nagel¹ & Jolyon P Mitchell²

¹Trudell Medical international, 725 Third Street, London, Ontario, N5V 5G4, Canada

²Jolyon Mitchell Inhaler Consulting Services Inc., 1154 St. Anthony Road, London, Ontario, N6H 2R1, Canada

Summary

The assessment through a challenge test of severity of reversible broncho-constrictive disease, such as asthma, is an important part of the diagnosis process as well as defining treatment therapy. Methacholine is frequently used as the inhaled challenge substance and is given by inhalation via a nebulizer for a fixed exposure time to the methacholine concentration. The challenge test involves progressively increasing the concentration of methacholine and measuring the forced expiratory volume in 1 s (FEV₁) after exposure at each concentration level. The test is terminated after the first instance at which FEV₁ decreases by more than 20% from the pre-test reference value. New recommendations from the European Respiratory Society (ERS) have recommended basing the result upon the delivered dose (µg) of methacholine causing a 20% fall in FEV₁, termed the provocative dose (PD₂₀), rather than the former metric of methacholine concentration (mg/mL), causing the same fall in FEV₁ (PC₂₀). Given the detail and complexity of the recent guidance, we follow a step-wise approach to explain each stage of the new bronchial challenge test, then illustrate how PD₂₀ is calculated and used to interpret the degree of airway hyper-responsiveness. Although any nebulizer with validated methacholine delivery data could be used to deliver the agent, we illustrate how to apply the methodology, based on the same breath-actuated nebulizer (AeroEclipse-II* BAN) as was used, through references, in the new guidance.

Introduction

The assessment through a challenge test of severity of reversible broncho-constrictive disease, such as asthma, is an important part of the diagnosis process as well as defining treatment therapy ^[1]. Methacholine is frequently used as the inhaled challenge substance, because the onset of symptoms upon inhalation of an appropriate concentration is rapid, and spontaneous recovery post-methacholine testing usually occurs within 45–60 min ^[2]. In practice, however, patients are usually given a bronchodilator at the end of testing to relieve challenge-induced bronchoconstriction more rapidly ^[2]. The bronchial challenge agent is given by inhalation via a nebulizer for a fixed exposure time to the concentration of methacholine. The provocation test involves progressively doubling the concentration of methacholine and measuring the forced expiratory volume in 1-s (FEV₁) after exposure at each concentration level. The test is terminated after the first instance at which FEV₁ decreases by more than 20% from the pre-test reference value.

New recommendations from an international European Respiratory Society task force have been published this year ^[3]. This technical standard, also endorsed by the American Thoracic Society, recommends basing the result of the bronchial challenge upon the delivered dose (mass expressed in µg) of methacholine causing a 20% fall in forced expiratory volume in 1 s (FEV₁). This is termed the provocative dose (PD₂₀), and replaces the former definition based on the provocative concentration of challenge agent resulting in a 20% reduction in FEV₁ (PC₂₀). This new end-point allows comparable results from either different aerosol delivery devices or protocols. Hence, the standard notes that any suitable nebuliser or dosimeter may be used, so long as the delivery characteristics are known ^[3].

It is recognized however that the change in approach to assess PD₂₀ rather than PC₂₀ has the potential to cause some confusion in how to execute the protocol in a practical manner. The purpose of the present interpretation is therefore to provide a simplified explanation with a practical, step-wise, example of how the test can be performed to meet the new standard.

Bronchial Challenge Testing – Drug Delivery System

The new standard allows for ‘any suitable nebulizer or dosimeter’ but does require characterization of the device output and particle size to be demonstrated. The example provided in this abstract uses data from the breath actuated nebulizer (AeroEclipse-II* BAN, Trudell Medical International, London, Canada) that is specifically referenced in the new 2017 standard, using independently reported tidal breathing data (both *in vitro* and *in vivo*). Such a breath-actuated device, that only delivers the medication when the patient inhales, has the additional benefit of affording minimal exposure of health care personnel to fugitive emissions ^[4], although a filter can be placed on the expiratory limb to eliminate such exposure altogether ^[3]. At least two independent clinical studies have recommended using this breath actuated nebulizer for methacholine challenge testing ^[4, 5].

How to perform the challenge test: Example calculation of PD₂₀

1) Prepare the Methacholine Solutions for Challenge test

The dilutions of methacholine concentrate can be prepared in the same way as with the previous 1999 guidance, prior to performing the challenge test and measurements of FEV₁. Table 1 shows an example of a schedule, based on the guidance in the new ERS document³.

Table 1 - Methacholine Concentrate Dilution Schedule in Which the Challenge Agent Concentration is Increased Four-Fold for Each Exposure

Label Mass of Concentrate (mg)	Start with:	Normal Saline Added to Effect Dilution (mL)	Obtain Diluted Concentration (C _{MC}) (mg/mL)	Code Letter to Provide Order of Dilution (see second column)
100	100 mg	+ 6.25	16.0	A
	3-mL of A	+ 9.0	4.0	B
	3-mL of B	+ 9.0	1.0	C
	3-mL of C	+ 9.0	0.25	D
	3-mL of D	+ 9.0	0.0625	E
	3-mL of E	+ 9.0	0.015625	F

2) Calculate the Delivered Doses at different Methacholine concentrations

In order to establish the delivered dose to the lungs (DD_{MC}) during a defined delivery duration, several key parameters regarding the nebulizer's output characteristics need to be known. For example, Appendix D of the new ERS standard³ provides the following information for the BAN:

- For 20 seconds of tidal breathing, the delivery rate (R_{MC}) of methacholine at the mouthpiece of the high output device (BAN) is 2.70 mg/min for a solution concentration (C_{MC}) of 16 mg/mL when operated from a 50-psi dry gas source.
- The fine droplet fraction (FDF), defined as those droplets less than 5 µm aerodynamic diameter, is reported from *in vitro* measurements of BAN-emitted droplets made by laser diffractometry as being 0.76³.

Hence the DD_{MC} for t(s) can be calculated as $DD_{MC} = R_{MC} \times FDF \times (t/60)$, and in the example provided for 20 seconds with the 16mg/mL concentration, DD_{MC} would therefore be:

- $2.70 \text{ mg/min} \times 0.76 \times 20/60 = 680 \text{ µg}$

This can further be generalized for any MC concentration using 20 seconds tidal breathing with the BAN as:

$$DD_{MC} = [C_{MC}/16 \text{ mg/ml}] \times 680 \text{ µg.}$$

3) Perform the bronchial challenge test

Once the calculations of DD_{MC} are completed for all the concentrations prepared as part of the test phase in Table 1, the measurement of FEV₁ can be conducted at increasing concentrations. Table 2 is an example of a bronchial challenge report taken from Appendix F of the ERS standard³. The DD_{MC} values in this case are based upon a 1 minute tidal breathing test duration as recommended in the standard.

The test begins with a 'Pre-Challenge' to confirm that the patient can perform acceptable and repeatable spirometry, and ensure they have sufficient airflow at baseline. Increasing amounts of DD_{MC} are delivered until such time as FEV₁ has fallen >20% from the reference (baseline) condition. In this particular example, in Table 2, the test was terminated after exposure to 127 µg (D₂) and the dose at the second to last exposure D₁ is 31.8 µg.

Table 2 - Example Bronchial Challenge Report

Time of exposure	Test Phase	DD _{MC} (µg) (1 min tidal breathing)	FEV ₁ (L)	FEV ₁ (% of reference)	Change in FEV ₁ (% pre-challenge value)
T ₀	Pre-challenge	N/A	3.10	N/A	N/A
T ₀ + 10 min	Post diluent (reference condition)	N/A	3.00	100	N/A
T ₀ + 15 min	0.015625 mg/mL	1.9	3.05	102	-2
T ₀ + 20 min	0.0625 mg/mL	7.65	2.94	98	2
T ₀ + 25 min	0.25 mg/mL	31.8	2.62	87	13
T ₀ + 30 min	1.0 mg/mL	127	2.16	72	28
T ₀ + 45 min	After bronchodilator administration	N/A	3.20	107	-7

4) Determination of PD₂₀

The PD₂₀ calculation is shown below and is illustrated using the example data from Table 2 where **R₁** and **R₂** are the percentage decreases in FEV₁ for **D₁** and **D₂**, respectively.

$$PD_{20} = \text{antilog} \left\{ \log D_1 + \frac{(\log D_2 - \log D_1)(20 - R_1)}{(R_2 - R_1)} \right\}$$

$$PD_{20} = \text{antilog} \left\{ 3.46 + \frac{(4.84 - 3.46)(20 - 13)}{(28 - 13)} \right\}$$

Consequently, from this particular example above, the bronchial responsiveness (PD₂₀) is determined as 61 µg.

5) Assessment of Airway Hyper-Responsiveness (AHR)

The PD₂₀ value can then be used to interpret the degree of AHR using values from the ERS document³ represented below in Table 3. Based on the given example, the patient would be considered to have Mild AHR.

Table 3 - Categorization of AHR to PD₂₀ of Methacholine

PD ₂₀ (µg)	Interpretation
>400	Normal
100–400	Borderline AHR
25–100	Mild AHR
6–25	Moderate AHR
<6	Marked AHR

Conclusions

The new ERS standard allows the use of a more appropriate PD₂₀ endpoint to assess airway hyper-responsiveness. The methacholine challenge test procedure, calculation and interpretation have been described in an attempt to provide a meaningful practical demonstration of how the new guideline could be put into practice clinically.

References

- ¹ Global Initiative for Asthma (GINA). Global strategy for asthma management and prevention. 2017 update. Available at: <http://ginasthma.org/2017-gina-report-global-strategy-for-asthma-management-and-prevention/> visited June 22 2017.
- ² Cockcroft DW, Swystun VA, Bhagat RG. Interaction of inhaled beta-2 agonist and inhaled corticosteroid on airway responsiveness to allergen and methacholine. *Am J Respir Crit Care Med* 1995; 152: pp1485–1489.
- ³ Coates AL, Leung K, Dell SD. Developing alternative delivery systems for methacholine challenge tests. *J Aerosol Med Pulmon Drug Deliv* 2014; 27: pp.66–70.
- ⁴ Dole SD, Bola SS, Foty RG, Marshall LC, Nelligan KA, Coates AL. Provocative dose of methacholine causing a 20% drop in FEV₁ should be used to interpret methacholine challenge tests with modern nebulizers. *Ann Am Thorac Soc*. 2015; 12(3): pp 357–363.
- ⁵ El-Gammal AI, Killian KJ, Scime TX, Beaudin S, Schlatman A, Cockcroft DW, Gauvreau GM. Comparison of the provocative concentration of methacholine causing a 20% fall in FEV₁ between the AeroEclipse II breath-actuated nebulizer and the Wright nebulizer in adult subjects with asthma. *Ann Am Thorac Soc*. 2015; 12(7): pp 1039–1043.