

The Methacholine Challenge Test for Reversible Airways Disease Assessment: A Practical Guide on How to Interpret New 2017 ERS Guidelines

Jason A. Suggett¹, Mark W. Nagel¹, Jolyon P. Mitchell²

¹Trudell Medical International, London, Canada. ²Jolyon Mitchell Inhaler Consulting Services Inc., London, Canada.

INTRODUCTION

- Methacholine is frequently used to assess the severity of reversible broncho-constrictive disease, such as asthma
- New ERS Technical Standards were published in 2017 recommending a shift from use of PC₂₀ to PD₂₀ for bronchial challenge testing
 - Allows comparable results from different aerosol delivery devices or protocols and allows use of any nebulizer where delivery characteristics are known
- The change from PC₂₀ to PD₂₀ has the potential to confuse how to execute the protocol in a practical manner
- The purpose of the present interpretation is therefore to provide a simplified explanation with a practical, step-wise, example of how the test can be performed to meet the new standard

BRONCHIAL CHALLENGE TESTING — DRUG DELIVERY SYSTEM

- The new standard allows for “any suitable nebulizer or dosimeter” with demonstrated device output and particle size characterization
- This poster will provide the data for the **AeroEclipse**® II Breath Actuated Nebulizer (BAN, Trudell Medical International, London, Canada) using independently reported tidal breathing data used in the new standard
- At least two independent clinical studies have recommended using this breath actuated nebulizer for methacholine challenge testing^{1,2}

HOW TO PERFORM THE CHALLENGE TEST: EXAMPLE CALCULATION OF PD₂₀

1 Prepare the Methacholine Solutions for Challenge Test

- Dilutions of methacholine concentrate can be prepared in the same way as with the previous 1999 guidance
- Table 1 shows an example of a schedule, based on the guidance in the new ERS document³

Table 1: Methacholine Concentrate Dilution Schedule in Which the Challenge Agent Concentration is Increased Four-Fold for Each Exposure

Label Mass of Concentrate (mg)	Start with	Normal Saline Added to Effect Dilution (mL)	Obtain Diluted Concentration (C _{MC}) (mg/mL)	Code Letter to Provide Order of Dilution (see second column)
100	100 mg	+6.25	16.0	A
	3 mL of A	+9.0	4.0	B
	3 mL of B	+9.0	1.0	C
	3 mL of C	+9.0	0.25	D
	3 mL of D	+9.0	0.0625	E
	3 mL of E	+9.0	0.015625	F

2 Calculate the Delivered Doses at different Methacholine Concentrations

- In order to establish the delivered dose to the lungs (DD_{MC}) during a defined delivery duration, the following need to be known
 - Delivery Rate of methacholine
 - Fine Droplet Fraction (droplets less than 5 µm aerodynamic diameter)
- Appendix D of the new ERS standard³ provides the following information for the BAN
 - For 20 seconds of tidal breathing, the delivery rate of methacholine (R_{MC}) at the mouthpiece of the high output device (BAN) is 2.70 mg/min for a solution concentration (C_{MC}) of 16 mg/mL when operated from a 50-psi dry gas source
 - The Fine Droplet Fraction (FDF), defined as those droplets less than 5 µm aerodynamic diameter, is reported from *in vitro* measurements of BAN-emitted droplets made by laser diffractometry as being 0.76³
- Hence the DD_{MC} for t(s) can be calculated as:

$$DD_{MC} = R_{MC} \times FDF \times (t/60)$$

- In the example provided for 20 seconds with the 16mg/mL concentration, DD_{MC} would therefore be: 2.70 mg/min X 0.76 X 20/60 = 680 µg
- This can further be generalized for any C_{MC} (concentration of methacholine), using 20 seconds tidal breathing with the BAN as:

$$DD_{MC} = [C_{MC} / 16 \text{ mg/ml}] \times 680 \mu\text{g}$$

3 Perform the bronchial challenge test

- Once the calculations of DD_{MC} are completed for all the concentrations prepared as part of the test phase in Table 1, the measurement of FEV₁ can be conducted at increasing concentrations.
- Table 2 is an example of a bronchial challenge report taken from Appendix F of the new standard³. The DD_{MC} values in this case are based upon a **1 minute tidal breathing test** duration as recommended in the standard.

Table 2: Example Bronchial Challenge Report

In this particular example, the test was terminated after exposure to **127 µg (D2)** and the dose at the second to last exposure **D1 is 31.8 µg**.

Time of exposure	Test Phase	DD _{MC} (µg) (1 min tidal breathing)	FEV ₁ (L)	FEV ₁ (% of reference)	Change in FEV ₁ (% pre-challenge value)
T ₀	Pre-challenge	N/A	3.10	N/A	N/A
T ₀ + 10 min	Post diluent (reference condition)	N/A	3.00	100	N/A
T ₀ + 15 min	0.015625 mg/mL	1.9	3.05	102	-2
T ₀ + 20 min	0.0625 mg/mL	7.65	2.94	98	2
T ₀ + 25 min	0.25 mg/mL	31.8	2.62	87	13
T ₀ + 30 min	1.0 mg/mL	127	2.16	72	28
T ₀ + 45 min	After bronchodilator administration	N/A	3.20	107	-7

4 Determination of PD₂₀

- The PD₂₀ calculation below illustrates use of the example data from Table 2 where R₁ and R₂ are the percentage decreases in FEV₁ for D₁ and D₂, respectively

$$PD_{20} = \text{antilog} \left\{ \log D_1 + \frac{(\log D_2 - \log D_1) (20 - R_1)}{(R_2 - R_1)} \right\}$$

$$PD_{20} = \text{antilog} \left\{ 3.46 + \frac{(4.84 - 3.46) (20 - 13)}{(28 - 13)} \right\}$$

- Consequently, from this particular example above, the **bronchial responsiveness (PD₂₀) is determined as 61 µg**

5 Assessment of Airway Hyper-Responsiveness (AHR)

- The PD₂₀ value can then be used to interpret the degree of AHR using values from the ERS document³ represented in Table 3

Table 3: Categorization of AHR to PD₂₀ of Methacholine

PD ₂₀ (µg)	Interpretation
>400	Normal
100–400	Borderline AHR
25–100	Mild AHR
6–25	Moderate AHR
<6	Marked AHR

- Based on the given example, the patient would be considered to have **Mild AHR**

CONCLUSIONS

- The new ERS standard allows the use of a more appropriate PD₂₀ endpoint to assess airway hyper-responsiveness
- The methacholine challenge test procedure, calculation and interpretation have been described in an attempt to provide a meaningful practical demonstration of how the new guideline could be put into practice clinically

REFERENCES

- Dole SD, Bola SS, Foly RG, Marshall LC, Nelligan KA, Coates AL. Provocative dose of methacholine causing a 20% drop in FEV₁ should be used to interpret methacholine challenge tests with modern nebulizers. *Ann Am Thorac Soc.* 2015; 12(3):357–363.
- El-Gammal AI, Killian KJ, Scime TX, Beaudin S, Schlatman A, Cockcroft DW, Gauvreau GM. Comparison of the provocative concentration of methacholine causing a 20% fall in FEV₁ between the **AeroEclipse**® II breath-actuated nebulizer and the Wright nebulizer in adult subjects with asthma. *Ann Am Thorac Soc.* 2015; 12(7): 1039–1043.
- Coates AL, et al. ERS technical standard on bronchial challenge testing: general considerations and performance of methacholine challenge tests. *Eur Respir J* 2017;49:1601526.

